UK – System 75 gas spring Wall Mounted EX


Fig. 1

For mounting in an EX working area the earthing wire on the socket pipe is very important.

Take care not to damage it during the mounting.

1. Fasten the wall bracket solidly to the wall or column. When mounting the wall bracket be aware of the considerable strain the operation of the arm will transmit to the bracket.
2. If the wall or column is made of a conductive material the connection to the wall bracket must be non-conductive.
3. Place the square-guide in the top of the bracket (1) and fasten it to the bracket.
4. When unfolding the extraction arm, make sure the socket pipe is turned in the direction of the arrow (2) and that the fingerscrews (I) are faced forward.
5. Place the socket pipe in the wall bracket (3) and fasten the square-guide with the included screws.
6. The extraction arm can be turned 190° in the wall bracket.


Fig. 2

1. Connect the main ventilation duct directly to the top of the wall bracket.
2. Unfold the extraction arm in the direction of the arrow (1). Push the extraction arm towards vertical position. This enables the gas spring to be mounted (2). Turn the gas spring towards the fitting on the socket pipe (3). Fasten it with the fingerscrew.
3. In an explosive working area a potential equalizing must be ensured for the ventilation system. It is therefore important to bond all the conductive parts of the system and to earth them. The extraction arm must be connected to the ventilation duct with the earthing wire pre-attached to the extraction arm. Alsident® System A/S recommends the earthing of the extraction arm to be carried out by personnel experienced in EX areas.
4. Place the connection pipe of the accessory in the joint. Turn the slit on the connection pipe towards the small fingerscrew on the joint. When screwed in the fingerscrew will hold the accessory in place in the joint.
5. Finally make sure the damper is open.
6. The extraction arm is now ready for use.

UK – System 75 gas spring Ceiling Mounted EX


Fig. 1

For mounting in an EX working area the earthing wire on the socket pipe is very important. Take care not to damage it during the mounting.

1. Mount the column directly under the ceiling.
2. The operation of the extraction arm will transmit a considerable strain to the column. The column must therefore be mounted on a solid ceiling.
3. When unfolding the extraction arm, make sure the socket pipe is turned in the direction of the arrow and the fingerscrews (I) are faced forward.
4. Place the socket pipe in the column and fasten the square-guide with the included screws.
5. The extraction arm can now be turned 360° in the column.


Fig. 2

1. For the side connection (Fig. 2a) the main ventilation duct (80mm dia.) can be connected directly to the outlet on the side of the column. For top connection (Fig. 2b) the main ventilation duct is put through the ceiling and directly into the column. Make sure to tighten the gap between the ceiling and duct with a filler.
2. Unfold the extraction arm in the direction of the arrow (1). Push the extraction arm towards vertical position. This enables the gas spring to be mounted (2). Turn the gas spring towards the fitting on the socket pipe (3). Fasten it with the fingerscrew.
3. In an explosive working area a potential equalizing must be ensured for the ventilation system. It is therefore important to bond all the conductive parts of the system and to earth them. The extraction arm must be connected to the ceiling column with the earthing wire pre-attached to the extraction arm. The installer of the ventilation duct must make sure to connect the ceiling column to the ventilations duct with an earthing wire. Alsident® System A/S recommends the earthing of the extraction arm to be carried out by personnel experienced in EX-areas.
4. Place the connection pipe of the accessory in the joint. Turn the small slit in the connection pipe towards the fingerscrew. When screwed in the fingerscrew will hold the accessory in place in the joint.
5. Finally make sure the damper is open.
6. The extraction arm is now ready for use.

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EX

UK: Approved for use in explosive areas and marked

Ex II 1 GD

according to standards:

EN 80079-36:2016 and EN 80079-37:2016, supporting ATEX directive 2014/34/EU.

For further information see Technical Information at: www.alsident.com.

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ESD

UK:

Approved according to standard IEC 61340-5-1:2016.

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